ENGLISH

PULP CANAL SEALER A ROOT CANAL LINING AND SEALER

FIELD OF APPLICATION

Pulp Canal Sealer is used for permanent obturation of the root canal space with the aid of obturating points in accordance with ISO 6876 for Dental Root Canal Sealing Materials.

SIDE EFFECTS

General

Zinc Oxide-Eugenol based sealers may produce periapical inflammation in the initial few days after placement. The inflammation subsides as healing progresses. A small number of patients may experience an allergic reaction to eugenol-containing materials, and caution is recommended.

DISPENSING AND APPLICATION

1. Procedure for measuring and dispensing:

Using either a mixing pad or a small glass slab, measure one level scoop of powder to mixing surface. Dispense one drop of eugenol liquid from the dropper, being careful to replace the caps of liquid and powder containers after use to keep humidity out.

2. Procedure for mixing:

Using a spatula, incorporate the powder into the liquid and mix in a very small area of about 1-2 cm. diameter. The mixing should be in a circular motion while pressing the spatula hard. The correct mix should be smooth and non-sticking.

3. Application:

Using Lentulo Spiral or appropriate Endodontic File, hand carry the mixed sealer into the canal. Coat the obturation point with the sealer using the spatula, then seat the point slowly.

4. Work and Set Time Specifications:

Pulp Canal Sealer WorkTime (Pad)	SetTime (Mouth)
70°-80°F (21°-27°C) Temperature	98.6°F (37°C)
40-60% Humidity	100%
> 45 Minutes Maximum	< 60 Minutes

5. Problem Solving

A. Mix is setting up too quickly:

Adding additional sealer liquid by moistening the tip of a spatula and incorporating it into the mix will lengthen the working time. Heat and humidity will accelerate the setting time of the mixed components. Any liquid dropper should be kept clean and dry, free of moisture contaminants. Check for age of materials, as working time decreases over time as materials pick up moisture while containers are open.

B. Mix is setting too slowly or not at all:

Check for contamination of the sealer liquid, powder, mixing surface, or liquid dropper.

STORE AND SHELF LIFE

Store packages in a cool and dry location, away from heating and high humidity sources.

PACKAGING AND PART NUMBERS

Pulp Canal Sealer:
Standard Package, Item# 4000-66
32g. Powder

WARRANTY INFORMATION/LIMITATION OF DOC'S BEST LIABILITY

Doc's Best technical advice, whether verbal or in writing, is designed to assist dentists in using Doc's Best's products. Such advice does not expand Doc's Best's limited warranty or relieve the dentist of testing Doc's Best products to determine their suitability for the intended uses and procedures. The dentist assumes all risks and liability for damages arising out of the improper use of Doc's Best's product.

In the event of a defect in material or workmanship, Doc's Best's liability is limited at Doc's Best's option, to replacement of the defective product or part thereof, or reimbursement of the actual cost of the defective product for a period covering one year from the date of purchase, with proof of purchase (original invoice) required. In order to take advantage of this limited warranty, the defective product must be returned to Cooley & Cooley. In no event shall Cooley & Cooley or Doc's Best be liable for any indirect, incidental, or consequential damages.

EXCEPT AS EXPRESSLY PROVIDED ABOVE, THERE ARE NO WARRANTIES BY DOC'S BEST EXPRESS OR IMPLIED, INCLUDING WARRANTIES WITH RESPECT TO DESCRIPTION, QUALITY OR FITNESS FOR A PARTICULAR PURPOSE.

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MATERIAL SAFETY DATA SHEET

SECTION I- PRODUCT IDENTIFICATION

PRODUCT NAME: Doc's Best Pulp Canal Sealer with Activated Copper

DISTRIBUTOR: Cooley & Cooley Ltd
Address: 8550 Westland West Blvd
Houston TX 77041-1214

FOR INFORMATION CALL: (281) 897-0009

INFOTRAC: 352-323-3500

DATE PREPARED: Feb. 5, 2004

SECTION II- HAZARDOUS INGREDIENTS/ IDENTIFICATION INFORMATION

HAZARDOUS COMPONENT	OSHA PEL	ACGIH TLV	%(OPTIONAL)
SODIUM BORATE	10 mg/m ³	5 mg/m ³	
ZINC OXIDE	10 mg/m ³	10 mg/m ³	
BARIUM SULFATE			
BISMUTH SUBCARBONATE			
HYDROGENATED ROSIN			
CUPROUS OXIDE			

SECTION III-PHYSICAL DATA

BOILING POINT: ND	SPECIFIC GRAVITY (H ₂ O = 1): ND
VAPOR PRESSURE (mmHg): ND	MELTING POINT: ND
VAPOR DENSITY (AIR=1): ND	EVAPORATION RATE: ND (BUTYL ACETATE=1)

SOLUBILITY IN WATER: NEGLIGIBLE

APPEARANCE AND ODOR: OFF WHITE SLIGHTLY GRANULAR POWDER

SECTION IV: FIRE AND EXPLOSION HAZARD DATA

FLASH POINT: ND	FLAMMABLE LIMIT	LEL: ND
		UEL: ND

EXTINGUISHING MEDIA: WATER SPRAY, DRY CHEMICAL, FOAM OR CARBON DIOXIDE.

SPECIAL FIRE FIGHTING PROCEDURES: FIREFIGHTERS SHOULD WEAR SELF-CONTAINED BREATHING APPARATUS.

UNUSUAL FIRE AND EXPLOSION HAZARDS: NONE

SECTION V – HEALTH HAZARD DATA

ROUTE (S) OF ENTRY –	INHALATION: POSSIBLE
	SKIN: POSSIBLE
	INGESTION: POSSIBLE

HEALTH HAZARD (ACUTE OR CHRONIC) ACUTE: MAY CAUSE MUCUS MEMBRANE IRRITATION. CHRONIC: REPEATED OR PROLONGED CONTACT MAY CAUSE DERMATITIS.

CARCINOGENICITY: ND NTP? IARC MONOGRAPHS? OSHA REGULATED? SIGNS OR SYMPTOMS OF EXPOSURE: SKIN SENSITIZATION MAY BE POSSIBLE IN SOME INDIVIDUALS.

MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE: ND

EMERGENCY AND FIRST AID PROCEDURES:

INHALATION: REMOVE TO FRESH AIR.

EYES: FLUSH WITH COPIOUS AMOUNTS OF WATER

SKIN: WASH WITH SOAP AND WATER

INGESTION: SEEK MEDICAL ATTENTION.

SECTION VI – REACTIVITY DATA

STABILITY: UNSTABLE: STABLE: X

CONDITIONS TO AVOID: ND

INCOMPATIBILITY (MATERIAL TO AVOID): ALUMINUM POWDER, PHOSPHORUS, AND ZIRCONIUM

HAZARDOUS DECOMPOSITION OR BYPRODUCTS: MAY INCLUDE TOXIC SODIUM OXIDE.

HAZARDOUS POLYMERIZATION: MAY OCCUR: WILL NOT OCCUR: X
CONDITIONS TO AVOID: ND

SECTION VII PRECAUTIONS FOR SAFE HANDLING AND USE

STEPS TO BE TAKEN IN CASE MATERIAL IS RELEASED OR SPILLED:

SWEEP UP ANY SPILLED MATERIAL.

WASTE DISPOSAL METHOD: DISPOSE OF IN ACCORDANCE WITH FEDERAL, STATE, AND LOCAL REGULATIONS

PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE: STORE MATERIAL AT ROOM TEMPERATURE AWAY FROM DIRECT SUNLIGHT.

OTHER PRECAUTIONS: ND

SECTION VIII- CONTROL MEASURES

RESPIRATORY PROTECTION (SPECIFY TYPE): ADEQUATE VENTILATION TO KEEP PARTICULATES BELOW EXPOSURE LIMITS:

VENTILATIONS – LOCAL EXHAUST: ND SPECIAL: ND
MECHANICAL (GENERAL): ND OTHER:ND

PROTECTIVE GLOVES: RECOMMENDED EYE PROTECTION: ND
OTHER PROTECTIVE CLOTHING OR EQUIPMENT: ND

SECTION X- ADDITIONAL INFORMATION

INFORMATION HEREIN GIVEN IN GOOD FAITH AS AUTHORITATIVE AND VALID. HOWEVER, NO WARRANTY, EXPRESSED OR IMPLIED CAN BE MADE, AND COOLEY & COOLEY LTD. ASSUMES NO LIABILITY RESULTING FROM ITS USE.