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

Bacterial filament

Diatom

ASM MicrobeLibrary.org © Donlan & Gibbons

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“Achieving Optimum Retention for Restorations”

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Cement	Copper Content	Bacterial Count Per CC. after 24 hours	Efficiency
Ames Black	97%	0	100%
Caulk's Black	25%	40,500	97.9%
Caulk's Red	8%	41,000	97.7%
Fleck's Red	7%	20,450	98.9%
Caulk's White	2%	81,000	95.9%
Smith's White	2%	81,000	95.9%

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- ❖ Evaluation of Cytotoxicity of the Mixture of Doc's Best Antibiofilm Cement Powder with Copalite Varnish Using the Agar Diffusion Method, Yiming Li, DDS., PhD., Wu Zhang, MD., Omari Onyango, DDS., MPH., Biocompatibility and Toxicology Research Laboratory, Loma Linda Univ. School of Dentistry, Sept. 2, 2005.

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Table 1
 Evaluation of Cytotoxicity of Mixture of Doc's Best Antibiofilm Cement Powder with Copalite Varnish Using the Agar Diffusion Method (24-hour Data)

Group ^a	Decolorization		Cell Lysis		Cell Response	Cytotoxicity
	cm ^b	Index	%	Index		
Negative Control	0	0	0	0	0/0	None
Mixture of Doc's Best with Copalite						
Fresh mixture	0	0	0	0	0/0	None
Set Material	0	0	0	0	0/0	None
Medium Extracts	0	0	0	0	0/0	None
Cottonseed Oil Extracts	0	0	0	0	0/0	None
Positive Control	1.45 ± 0.00	5	64.0 ± 0.04	4.0	5/4	Severe

^a Mixture of Doc's Best with Copalite: the mixture of Doc's Best Antibiofilm Cement Powder with Copalite Varnish. Negative Control: filter disk with 50 µL sterile culture media without serum; Positive Control: filter disk with 50 µL phenol.

^b N=4. The distance from the sample (cm) = (Diameter of the Decolorization Zone - Diameter of the sample) / 2. The value of 1.45 cm indicates a decolorization of entire culture well (3.5 cm in diameter): 1.45 cm x 2 + 0.6 cm (diameter of the sample). Decolorization Index is 1 if the Decolorization Zone is limited to the area under the sample (Appendix A).

Table 2
 Evaluation of Cytotoxicity of Mixture of Doc's Best Antibiofilm Cement
 Powder with Copalite Varnish Using the Agar Diffusion Method
 (48-hour Data)

Group ^a	Decolorization		Cell Lysis		Cell Response	Cytotoxicity
	cm ^b	Index	%	Index		
Negative Control	0	0	0	0	0/0	None
Mixture of Doc's Best with Copalite						
Fresh mixture	0	0	0	0	0/0	None
Set Material	0	0	0	0	0/0	None
Medium Extracts	0	0	0	0	0/0	None
Cottonseed Oil Extracts	0	0	0	0	0/0	None
Positive Control	1.45 ± 0.00	5	90.0 ± 0.04	5	5/5	Severe

^a Mixture of Doc's Best with Copalite: the mixture of Doc's Best Antibiofilm Cement Powder with Copalite Varnish. Negative Control: filter disk with 50 µL sterile culture media without serum; Positive Control: filter disk with 50 µL phenol.

^b N=4. The distance from the sample (cm) = (Diameter of the Decolorization Zone - Diameter of the sample) / 2. The value of 1.45 cm indicates a decolorization of entire culture well (3.5 cm in diameter): 1.45 cm x 2 + 0.6 cm (diameter of the sample). Decolorization Index is 1 if the Decolorization Zone is limited to the area under the sample (Appendix A).

Appendix A

Evaluation Criteria for Agar Diffusion Test

A. Zone Index

0	No detectable zone around or under sample
1	Zone limited to area under sample
2	Zone <0.5 cm from sample
3	Zone <1.0 cm from sample
4	Zone >1.0 cm from sample but not entire plate
5	Zone involves entire plate

B. Lysis Index

0	No observable lysis
1	<20%
2	20 - 40% of zone lysed
3	41 - 59% of zone lysed
4	60 - 80% of zone lysed
5	>80% of zone lysed

C. Procedures

1. The area of discoloration will be evaluated and measured for each plate.
2. The extent of lysis within a zone will be estimated using an inverted microscope.

D. Response Index

Cell Response = Decolorization Index / Lysis Index

E. Interpretation of the results

Cell Response	Score	Interpretation
0/0 - 0.5/0.5	0	Nontoxic
1/1 - 1.5/1.5	1	Mildly toxic
2/2 - 3/3	2	Moderately toxic
4/4 - 5/5	3	Severely toxic

Evaluating the Anti-Plaque Capabilities of a Copper Containing Prophylaxis Paste, Moore, Feldman, Abbott and Associates, L.S.U. School of Dentistry, 2004

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- ❖ **Red Copper Cement**
- ❖ **White Copper Cement**
- ❖ **Red Pulp Canal Sealer**
- ❖ **White Pulp Canal Sealer**

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Table 1. Test tube designation, inoculum, media and average colonization.

Tube Numbers	Coupon Type	Innoculum	Media	Average Colonization (CFU/cm ²)
1-3	GLS	<i>S. mutans</i>	BHI	0
4-6	GLS	<i>Lactobacillus</i>	MRS	684
7	GLS	none	BHI	0
8	GLS	none	MRS	0
9-11	RCC	<i>S. mutans</i>	BHI	0
12-14	RCC	<i>Lactobacillus</i>	MRS	0
15	RCC	none	BHI	0
16	RCC	none	MRS	0
17-19	WC	<i>S. mutans</i>	BHI	0
20-22	WC	<i>Lactobacillus</i>	MRS	0
23	WC	none	BHI	0
24	WC	none	MRS	0
25-27	WC+	<i>S. mutans</i>	BHI	0
28-30	WC+	<i>Lactobacillus</i>	MRS	0
31	WC+	none	BHI	0
32	WC+	none	MRS	0
33-35	PC	<i>S. mutans</i>	BHI	0
36-38	PC	<i>Lactobacillus</i>	MRS	210
39	PC	none	BHI	0
40	PC	none	MRS	0

Evaluation of Mutagenic Potential of the Mixture of Doc's Best White Copper Activated Copper Cement with Copalite Varnish Using the Ames *Salmonella*/Microsome Mutagenicity Test

Study conducted at Loma Linda University School of Dentistry, Biocompatibility and Toxicology Research Laboratory, Submitted December 22, 2005.

Researchers: Yiming Li, D.D.S., Phd., Wu Zhang, M.D., Minling Zheng, D.D.S., M.S.

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Table 1
 Evaluation of Mutagenic Potential of the Mixture of Doc's Best White Copper Activated
 Copper Cement with Copalite Vanish Using the Ames *Salmonella* Mutagenicity Test
 (Without S9 Activation)

Test Agent ^a	Dose (Dilution)	Revertants/Plate			
		TA97a	TA98	TA100	TA102
Negative Control	—	123 ± 4	18 ± 4	169 ± 12	195 ± 30
Doc's Best/Copalite Extracts	1:9	98 ± 6	14 ± 4	169 ± 6	191 ± 1
	1:4	92 ± 0	19 ± 1	159 ± 1	218 ± 16
	1:2	88 ± 11	12 ± 4	187 ± 29	201 ± 24
	As is	103 ± 13	11 ± 1	180 ± 1	211 ± 6
Positive Control	—	1004 ± 130	2028 ± 504	1172 ± 127	1536 ± 91

^a Doc's Best/Copalite: mixture of Doc's Best White Copper Cement with Copalite Vanish; Negative Control: sterile saline; Positive controls: TA97a and TA98, 4-nitro-o-phenylenediamine (20 µg/plate); TA100, Sodium azide (5 µg/plate); TA102, t-butyl hydroperoxide (100 µg/plate).


^b Mean ± standard deviation (N=2). Values within brackets are not significantly different as determined using the Student-Newman-Keuls method.

Table 2
 Evaluation of Mutagenic Potential of the Mixture of Doc's Best White Copper Activated
 Copper Cement with Copalite Vanish Using the Ames *Salmonella* Mutagenicity Test
 (With S9 Activation)

Test Agent ^a	Dose (Dilution)	Revertants/Plate			
		TA97a	TA98	TA100	TA102
Negative Control	—	257 ± 42 ^b	21 ± 4	143 ± 16	353 ± 24
Doc's Best/Copalite Extracts	1:9	195 ± 47	17 ± 3	150 ± 9	332 ± 2
	1:4	251 ± 35	19 ± 4	141 ± 1	290 ± 3
	1:2	218 ± 0	19 ± 5	163 ± 1	358 ± 28
	As is	227 ± 23	17 ± 6	169 ± 14	360 ± 54
Positive Control	—	1276 ± 96	3304 ± 238	1944 ± 11	1945 ± 69

^a Doc's Best/Copalite: mixture of Doc's Best White Copper Cement with Copalite Vanish; Negative Control: sterile saline; Positive control: TA97a, TA98 and TA100, 2-aminofluorene (15 µg/plate); TA102, 1,8-dihydroxyanthraquinone (50 µg/plate).

^b Mean ± standard deviation (N=2). Values within brackets are not significantly different as determined using the Student-Newman-Keuls method.

Applications cliniques

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
Tooth # 18- Through and through decay
(x-ray taken 2-11-99)

Six ans

Tooth # 18- Two years later


Cas clinique

**Tooth # 18- Composite filling removed
after two years displaying red copper cement
indirect pulp cap**

**TOOTH # 18- Secondary Dentin Formation
after two years**

Tooth # 18- Vital, no root calcification
(both x-rays taken on 2-21-05)

Après six ans

Tooth # 18- Vital, no root calcification
(both x-rays taken on 2-21-05)

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
Deux ans


Tooth # 13- Decay through to pulp chamber
(x-ray taken in August 2003)


Tooth #13- Vital with secondary dentin formation
Doc's Best Red Copper base (x-ray taken, August 2005)

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
Applications Endodontiques

ENDODONTIC TECHNIQUES

Placement of endodontic posts and replacement of failing posts

Step 1

If root canal filling is not sound, we recommend retreating the tooth endodontically using DOC'S BEST™ Pulp Canal Sealer and gutta percha.

Step 2


A post space is prepared using appropriate tools to remove either a failing previously placed post or by removing gutta percha.



NOTE: DOC'S BEST White Cement is outstanding to treat internal root resorption and problems with re-decay. Phosphoric acid cements are not indicated for perforations of the pulpal floor of teeth.


All Doc's Best™ cements are biofilm resistant, however, phosphate cements, ideal for indirect pulp caps and cementations, are inappropriate for repairing perforations in the pulp chambers and roots. A thick mix of Doc's Best™ Pulp Canal Sealer may be used to repair perforations.

Step 5

Using air and a brush, saturate the post space with a slurry mixture of 4-8 drops of COPALITE® Varnish and ¼ scoop of Copper Cement Powder. Let dry for one minute to insure maximum penetration of dentinal tubules. Then dry with air syringe. COPALITE® Varnish itself is very antimicrobial. COPALITE® has the unique penetrating ability needed to deliver the activated copper ions into the dentin tubules offering unparalleled protection from Biofilm formation.

A decorative header at the top of the slide. It features a central white globe with a blue and white pattern, set against a background of a stylized landscape with green hills and a blue sky. The globe is flanked by two identical rectangular panels, each containing the same stylized landscape. The entire header is set against a dark red background.

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
TREATING GERIATRIC AND SPECIAL NEEDS PATIENTS

Buccal


Interproximal

Lingual


Step 1

The preparation is thoroughly saturated with several layers of COPALITE® Varnish (4-8 drops) and DOC'S BEST™ cement powder (1/4 scoop) in a slurry mixture. Decay removal in badly neglected teeth is often relative. Our goal must be to contain decay long enough to let the patient outlive their teeth.

Step 2

Depending on the size of the decayed area, several methods of cement application may be used. A plastic matrix may be held in position with either fingers or a matrix holder. The same is true for metal bands. The worst case scenario would be subgingival decay in a hemorrhagic area which circumscribes a 360° area of the tooth.

Two holes are drilled in metal matrix band and/or perforated in a plastic matrix band, one on the mesial-labial and the other on the distal-labial aspect.

DOC'S BEST™ Cement is mixed to a thick but still flowable consistency and injected with a CENTRIX® Tube and Plug into either the mesial-labial hole or the distal-labial hole until the cement is freely expressed from the opposite hole.

At that point finger coverage of the first hole of application is used and the syringe is moved to the opposite hole which is further injected until the cavity is filled to completion. The area is covered with cotton and allowed to set for five minutes. A finishing bur or diamond can be used for final contouring.

Another method of inserting the cement, is to dip an amalgam plugger in the cement powder, then fill with the thick mix of cement. Pack in the cement like amalgam.

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